



MediSafe – Stay away and defeat diseases

2022 - 143

Our Team



Supervisor

Mr. Ravi Supunya



Co- Supervisor

Mr. Samantha Rajapaksha



External Supervisor

Dr. Shashika
Liyanage



IT18077698

Thennakoon
T.M.B.C.K



IT19015422

Perera B.A.A.W.S



IT19015040

Rasuni Wageesha
H.A



IT19011608

Senanayaka
S.A.M.A.B.M

Introduction

- ❑ What is MediSafe ?
- ❑ Why are we doing this research?
- ❑ What is the purpose of our research?



Research problem

- There are some diseases that have arisen at present. (Heart attack, Pneumonia, Wheezing, Dengue, Covid'19)
- High cost for diagnosis.
- Informal lifestyle and busyness.
- Don't have enough idea about current situation of the country.

Background

Table 22. Leading Causes of Hospital Deaths, 2010 - 2019

Disease and ICD (10 th Revision) Code	2019		2018		2017		2016		2015		2014		2013		2012		2011 ¹		2010 ²	
	Rank	%	Rank	%	Rank	%	Rank	%	Rank	%	Rank	%	Rank	%	Rank	%	Rank	%	Rank	%
Ischaemic heart disease I20 - I25	1	15.1	1	15.0	1	14.2	1	14.1	1	14.2	1	14.8	1	14.7	1	14.4	1	13.4	1	12.8
Zoonotic and other bacterial diseases A20 - A49	2	12.1	3	10.9	2	11.5	3	11.6	3	9.7	3	9.1	6	7.9	6	7.1	6	6.7	6	6.6
Neoplasms ¹ C00 - D48	3	11.7	2	11.7	3	10.5	2	12.0	2	11.0	2	11.7	2	11.2	2	11.6	2	11.8	2	11.1
Diseases of the respiratory system excluding diseases of upper respiratory tract, pneumonia and influenza I20 - J22, J40 - J98	4	10.7	4	9.9	4	9.8	5	8.3	4	9.2	6	8.0	5	7.9	5	7.2	5	6.9	5	7.0
Pneumonia J12 - J18	5	8.0	7	7.8	6	8.2	7	6.4	7	7.5	7	6.6	8	6.1	8	5.7	9	5.2	9	5.2
Pulmonary heart disease and diseases of the pulmonary circulation I26 - I51	6	7.6	6	7.9	5	8.5	4	8.7	5	8.3	4	8.6	4	8.4	3	9.0	4	8.7	3	8.7
Cerebrovascular disease I60 - I69	7	7.6	5	8.0	7	7.7	6	8.2	6	8.2	5	8.4	3	8.6	4	8.7	3	8.7	4	8.7
Diseases of the urinary system N00 - N39	8	5.8	8	5.8	8	5.9	8	6.3	8	6.2	8	6.3	7	6.2	7	6.3	7	5.7	8	5.7
Diseases of the gastro-intestinal tract K20 - K92	9	5.0	9	5.1	9	5.1	9	5.5	9	5.3	9	5.7	9	5.7	9	5.4	8	5.4	7	6.2
Traumatic injuries S00 - T19, W54	10	3.6	10	3.9	10	3.8	10	3.9	10	3.8	10	3.5	11	3.3	11	3.7	11	3.6	11	3.7
Disease of the nervous system G00 - G98	11	1.3	13	1.4	14	1.4	14	1.4	17	1.3	16	1.4	15	1.4	16	1.5	19	1.4	18	1.6
Symptoms, signs and abnormal clinical and labo R00 - R99	12	1.3	11	1.5	12	1.5	12	1.6	11	2.3	11	3.2	10	4.8	10	4.5	10	4.1	10	5.0
Diabetes mellitus E10 - E14	13	1.3	12	1.4	11	1.7	11	1.8	13	1.6	13	1.6	13	1.6	14	1.7	14	1.9	16	1.7

¹ Includes deaths reported from the Cancer Hospital (not analysed by site and type of neoplasm)

² Excludes Mullativu District

Source : Medical Statistics Unit



Source : Medical Statistics Unit Ministry of Health
http://www.health.gov.lk/moh_final/english/public/elfinder/files/publications/AHB/AHS%202019.pdf

Overall solutions – 100%

- Developed an Arduino-based device that detects certain types of symptoms to diagnose certain heart and lung related diseases.
- Use some machine learning based techniques to identify diseases and clarify it.
- Show diseases spread rate to the user.
- Developing a web application to facilitate patient usage.



Research Objectives




Implement a device to get parameters of the patient and identify Covid'19. (Possibility as a percentage)

Disease level wise identification and provide suggestions/ recommendations to reduce the risk level.

Identify the exact lung disease among other lung diseases.

Identify the three major diseases spread rate in Sri Lanka.

System Diagram



Focusing areas




Hospital / Medical centers

Schools



Offices

Crowded areas





IT18077698

Thennakoon T.M.B.C.K

Specialization | Information Technology




Table of Contents

01

Research question

02

System diagram

03

Progress

04

Technologies

05

Requirements

06

Risk mitigation

07

Completion of the
project




Research question

- ❑ Identify all measurements using single device with few minutes.
- ❑ Simple and user-friendly web application and mobile application.




- ❑ Provide probability to infect Covid 19 & Give what are the necessary actions need to get by patient.
- ❑ Get necessary inputs and Generate healthy recommendations to day-to-day life.

System diagram




Progress




Latest technologies in MediSafe

MediSafe device



ESP8266
Max30105sens
Flow meter
Firebase

Dataset collection and data training



Kaggle
Regression
Python / React

IDE



Jupyter Notebook
Arduino IDE
Visual studio code
Anaconda Prompt

Requirements

Functional

- Interoperability
- Accuracy
- Compliance

Non – functional

- Maintainability
- Manageability
- Usability
- Integrity




Risk mitigation



- This is related to medical industry research therefore the accuracy should be compulsory.
- Need to compare the actual output with a recommendation of a doctor
- Need to confirm that the output should provide the correct according to the currently available devices.

Completion of the project



Medisafe
device

```
final2 | Arduino 1.8.19 (Windows Store 1.8.57.0)
File Edit Sketch Tools Help

final2 $  Firebase  flow_meater  max30105  so2  temp

16 int count = 0;
17 int temp_bps;
18 //-----Firebase-----
19
20 #include <ArduinoJson.h>
21 #include "FirebaseESP8266.h"
22 #include <ESP8266WiFi.h>
23 // Set these to run example.
24 #define FIREBASE_HOST "medisafe-research-default-rtdb.firebaseio.com/unit_1"
25 #define FIREBASE_AUTH "qjnABtFp7TrCsENApoxBQSe12lkghAo10PwrBB5"
26 #define WIFI_SSID "supun"
27 #define WIFI_PASSWORD "supun111191"
28 FirebaseData firebaseData;
29
30 #define SENSOR D4
31 long currentMillis = 0;
32 long previousMillis = 0;
33 int interval = 1000;
34 //boolean ledState = LOW;
35 float calibrationFactor = 4.5;
36 volatile byte pulseCount;
37 byte pulseSec = 0;
38 float flowRate;
39 unsigned int flowMilliLitres;
40 unsigned long totalMilliLitres;
41 // -----led
42
43 #define REDLED D5
44 // -----led
```

Arduino IDE

Completion of the project



```
C:\WINDOWS\system32\cmd.exe


(covid) C:\Users\user>cd C:\Users\user\Desktop\24-04-2022\covid

(covid) C:\Users\user\Desktop\24-04-2022\covid>C:


(covid) C:\Users\user\Desktop\24-04-2022\covid>python Runcovid.py
type oxygen level : 90
type your pulse : 96
type your Temperature : 90
confidence : 100.0 %
The probability of having a covid infection is 35.360000000000004%
Traceback (most recent call last):
  File "Runcovid.py", line 1, in <module>
    from covid import predictc
ImportError: cannot import name 'predictc' from 'covid' (C:\Users\user\Desktop\24-04-2022\covid\covid.py)

(covid) C:\Users\user\Desktop\24-04-2022\covid>
```

Output of
prediction



Completion of the project

 Medisafe Health System

[Home](#) [Contact Us](#)

Covid

Enter Oxygen Level

88

Enter Pulse

81

Enter Temperature

104



Risk of Covid Infection

82.96000000000004

View Suggestions

mediSafe.research@gmail.com

Copy Rights @2022



IT19015422

Perera B.A.A.W.S

Specialization | Information Technology




Table of Contents

01

Research question

02

System diagram

03

Progress

04

Technologies

05


Requirements

06

Risk mitigation

07

Completion of the
project



Research question

- How to identify the people who are suffering in such lung and heart diseases (level wise)


- How to check current situation in cost effectively

- How to provide easily recommendations via web application to the user


- What are the solutions we can give due to shortage of medicines



System diagram



Progress



Latest technologies in MediSafe

DataSet collection

kaggle



- <https://data.world/informatics-edu/heart-disease-prediction>
- <https://www.kaggle.com/datasets/johnsmith88/heart-disease-dataset>

Model training and Language



- K – nearest neighbor algorithm
- Libraries – pandas, sklearn, joblib, numpy
- Python / React

IDE

ANACONDA



- Vs Code
- Jupyter notebook
- Anaconda prompt

Requirements

Functional

- Interoperability
- Authentication.
- Report generate
- User friendly

Non - functional

- Quality
- Durability
- Security
- Privacy



Risk mitigation



- ✓ Entering current situation features difficult to known by person . So that those features will get from the implemented device.
- ✓ Adults are not well fluent in new technologies.
- ✓ Validity of the disease level will depend on the user inputs.

Completion of the project

Frontend

1

```
</div>
<LayoutHeader />
{loading && <$Spin />}
<$Row style={{ marginTop: "100px", marginLeft: "5%" }}>
  <h1>Pneumonia</h1>
</Row>
<$Row style={{ marginRight: "5%" }}>
  <$Col xl={12} sm={12}>
    <Row className="jus-con-cen row-items">
      <$Col xl={5} sm={12}>
        <p>Enter Age</p>
      </$Col>
      <$Col xl={10}>
        <Input
          name="age"
          handleChange={this.onHandleChange}
          value={form.age}
        />
      </$Col>
    </Row>
    <Row className="jus-con-cen row-items">
      <$Col xl={5} sm={12}>
        <p>Enter Gender</p>
      </$Col>
      <$Col xl={10}>
        <Radio.Group
          value={form.Gender}
          onChange={(e) => {
            this.onHandleChange("Gender", e.target.value);
          }}
        />
      </$Col>
    </Row>
  </$Col>
</Row>
```

Completion of the project

API controller implementation

1

```
@app.route('/risk', methods=['GET', 'POST'])
def predictR():
    data = {}
    post_data = request.json


    age = str(post_data['age'])
    Gender = str(post_data['Gender'])
    Cholesterol = str(post_data['Cholesterol'])
    Pulse = str(post_data['Pulse'])
    Smoke = str(post_data['Smoke'])
    Alcohol = str(post_data['Alcohol'])
    Shortness_of_breath = str(post_data['Shortness_of_breath'])
    Anxiety = str(post_data['Anxiety'])

    y_predictH, y_predictP, y_predictW = get_risk_level(
        age, Gender, Cholesterol, Pulse, Smoke, Alcohol, Shortness_of_breath, Anxiety)

    data['prediction_heart'] = y_predictH[0]
    data['prediction_wheeze'] = y_predictP[0]
    data['prediction_pneumonia'] = y_predictW[0]

    return jsonify(data)
```

Completion of the project

 Medisafe Health System

[Home](#) [Contact Us](#)

Pneumonia

Enter Age

34

Enter Gender

☐ Male ☒ Female

Enter Pulse

90

Smoke

☐ Yes ☒ No

Anxiety

☐ Yes ☒ No

Alcohol Usage

☒ Yes ☐ No

Shortness of Breath



☒ Yes ☐ No

Risk of pneumonia

Low

View Suggestions

1. Stay hydrated. Drink plenty of fluids, especially water, to help loosen mucus in your lungs.
2. Take your medicine as prescribed. Take the entire course of any medications your doctor prescribed for you. If you stop taking medication too soon, your lungs may continue to harbor bacteria that can multiply and cause your pneumonia to recur.
3. Check oxygen saturations and provide supplemental oxygen if saturations are <90%
4. Stop smoking – Smoking increases your risk for pneumonia and other health conditions. If you are a smoker, consider stopping.
5. Get lots of rest- Rest will help your body fight the infection.
6. Drink plenty of fluids. Fluids will keep you hydrated. They can help loosen the mucus in your lungs. Try water, warm tea, and clear soups.
7. Use a cool-mist humidifier or take a warm bath. This will help clear your

mediSafe.research@gmail.com Copy Rights @2022  



IT19015040

Rasuni Wageesha H.A

Specialization | Information Technology




Table of Contents

01

Research question

02

System diagram

03

Progress

04

Technologies

05

Requirements

06

Risk mitigation

07

Completion of the
project





Research question




There are many different types of lung diseases and diagnosing one might be difficult.

System diagram





Technologies in MediSafe



Requirements

Functional

- Upload the lung image to the system.

Non – functional

- Performance
- Availability
- Reliability



Risk mitigation



Validity of the clear images.

Completion of the project

API

```
1 @app.route('/lung', methods=['GET', 'POST'])
2 def predict():
3     data = {}
4
5     filestr = request.files['file'].read()
6     img = imread(filestr)
7
8     prediction_inst = []
9     prediction_conf = []
10    list_of_cf = []
11
12    outputs = model_process_img(img)
13    for item in outputs['predictions']:
14        list_of_cf.append(item['confidence'])
15
16    for item in outputs['predictions']:
17        if item['confidence'] == max(list_of_cf):
18            print(item['label'], max(list_of_cf) * 100)
19            conf = max(list_of_cf) * 100
20            prediction_inst.append(item['label'])
21            prediction_conf.append(conf)
22
23    temp_val = prediction_inst[0]
24    temp_conf = prediction_conf[0]
25    print(temp_val)
26
27    prediction_inst.clear()
28    prediction_conf.clear()
29    list_of_cf.clear()
30    data['detection'] = temp_val
31    data['detectionScore'] = temp_conf
32
33    return jsonify(data)
```

Completion of the project

The screenshot displays the MediSafe Health System web application. The browser address bar shows `localhost:3000/disease9`. The application header includes the MediSafe logo and navigation links for Home and Contact Us. The main content area is titled "Upload your X-ray to Identify disease". On the left, there is a file upload area with a dashed border and a blue folder icon. Below the upload area, a file named "Lung_Opacity-3.png" is shown with a small thumbnail. On the right, the results are displayed: "Disease" is "Lung Opacity" and "Confidence" is "99.99128580093384". A "View Suggestions" button is present. Below the results, a list of recommendations is shown: "Don't Smoke.", "Breathing Exercise Regularly.", "Get Regular Check-ups.", and "Avoid Exposure to Indoor Pollutants That Can Damage Your Lungs." The footer contains the email `mediSafe.research@gmail.com`, the copyright notice "Copy Rights @2022", and social media icons for Facebook and YouTube.

MediSafe Health System

Home Contact Us

Upload your X-ray to Identify disease

Click or drag file to this area to upload
Please upload X-ray images

Lung_Opacity-3.png

Disease: Lung Opacity
Confidence: 99.99128580093384

View Suggestions

- Don't Smoke.
- Breathing Exercise Regularly.
- Get Regular Check-ups.
- Avoid Exposure to Indoor Pollutants That Can Damage Your Lungs.

mediSafe.research@gmail.com

Copy Rights @2022



IT19011608

Senanayaka S.A.M.A.B.M

Specialization | Information Technology




Table of Contents

01

Research question

02

System diagram

03

Progress

04

Technologies

05

Requirements

06

Risk mitigation

07

Completion of the
project




Research question




Identify the disease count on the Sri Lanka.

Target Domain


System diagram



Progress



Technologies in MediSafe



- 1) <https://www.healthdata.org/sri-lanka>
- 2) <https://www.kaggle.com/>

Long short term-memory algorithms

Python / React

Jupyter notebook and vs code

Requirements

Functional

- Identify the spread rate count
- Display the data healthcare dashboard

Non – functional

- Accuracy
- Availability



Risk mitigation



Find relevant and accurate data.

Completion of the project

API

```
1 @app.route('/spread', methods=['GET', 'POST'])
2 def predictS():
3     data = {}
4     post_data = request.json
5
6     sickness = str(post_data['sickness'])
7     city = str(post_data['city'])
8     Date = post_data['Date']
9
10    print(sickness, city, Date)
11
12    responseArray = []
13    for x in Date:
14        response = get_prediction(sickness, city, str(x).split('T')[0])
15        responseArray.append(response)
16
17    data['details'] = responseArray
18    return jsonify(data)
```

Completion of the project

Frontend-submit detail


```
1 submit = async () => {
2   const form = { ...this.state.form };
3   const sickness = form.sickness;
4   const city = form.city;
5   let getDaysArray = this.getDaysArray(
6     moment(form.Date[0]).format("YYYY-MM-DD"),
7     moment(form.Date[1]).format("YYYY-MM-DD")
8   );
9
10  console.log(getDaysArray, "getDaysArray");
11
12  const data = {
13    sickness: sickness,
14    city: city,
15    Date: getDaysArray,
16  };
17  this.setState({ loading: true, getDaysArray: getDaysArray });
18  try {
19    await fetch("/spread", {
20      method: "POST",
21      headers: {
22        "Content-Type": "application/json",
23      },
24      body: JSON.stringify(data),
25    })
26    .then((response) => response.json())
27    .then((response) => {
28      var spread = response["details"];
29      this.setState({
30        spread: spread,
31      });
32    });
33    this.setState({ loading: false });
34    let graph = this.getGraph(this.state.spread, this.state.getDaysArray);
35    this.setState({
36      graph: graph,
37    });
38  } catch (error) {
39    this.setState({ loading: false });
40  }
41  this.setState({ loading: false });
42  };
```

Completion of the project

Frontend-graph

```
1 {graph && (  
2     <Line  
3         data={graph}  
4         xField="Date"  
5         yField="scales"  
6         autoFit={false}  
7         xAxis={{ range: [0, 1], tickCount: 5 }}  
8         point={{ size: 5, shape: "diamond" }}  
9         label={{  
10             style: {  
11                 fill: "#aaa",  
12             },  
13         }}  
14     />  
15 )}
```

Completion of the project



Commercialization



Commercialization



01

We discussed with some surgeons, and they allow us to test that implemented system (100%) in their premises. (Dispensary, Hospitals)

02

After completing (100%) our project we are hoping to introduce this to some clinical centers.


03

Advertise our mediSafe product via social media with its benefits.


The background features two large, abstract, organic shapes in a medium blue color. One shape is on the left side, curving upwards and towards the center. The other is on the top right, curving downwards and towards the center. They frame the central text.

Demonstration

Device implementation



Web app implementation





Thank You

Do you have any questions?